

Philips Sonicare i InterCare Standard sonic toothbrush heads

2-pack

Standard size Click-on BrushSync mode pairing

HX9002/10



Deeper reach between teeth. Superior cleaning.*

Clean between teeth and in hard-to-clean areas

With advanced tufting technology and extra-long bristles, i InterCare helps to remove up to 7x more plaque between teeth than a manual toothbrush and improve gum health in just two weeks.

For exceptional plaque removal

• Remove up to 7x more plaque with our InterCare brush head*

For exceptional gum care

• Improves gum health in just 2 weeks

Innovative technology

- Automatically selects optimal mode for ultimate results**
- Creates microbubbles that reach deep between teeth

Designed to optimize performance

- Always know when to replace. Always get an effective clean.
- Click-on design for simple brush head placement
- · Tested to meet your oral health needs



Highlights

Remove up to 7x more plaque



Click on our InterCare brush head to improve gum health in just 2 weeks. Extra-long bristles help remove more plaque in hard-to-reach spots and between teeth for healthy gums.

Extra-long bristles



The InterCare brush head features extra-long, densely-packed, high-quality bristles to target hidden plaque caught deep between your teeth. It is clinically proven to reduce gingival bleeding and inflammation in just 2 weeks.

BrushSync™ mode-pairing



You'll always get the best possible clean with our BrushSync™ mode-pairing feature**. The

InterCare brush head syncs with your BrushSync™ enabled Philips Sonicare toothbrush handle**, selecting the optimal brushing mode and intensity level for exceptional cleaning. All you need to do is start brushing.

Brush head replacements



Brush heads become less effective after 3 months of use, but with BrushSync™ you'll be reminded before this happens. Your smart toothbrush will track how often and how hard you brush, and will notify you when it's time for a replacement. Don't have a smart Philips Sonicare toothbrush? Just monitor your blue replacement bristles, and when they turn white you'll know it's time for a fresh brush head.

Easy to click-on



Your InterCare brush head fits perfectly with any Philips Sonicare toothbrush handle, except PowerUp Battery and Essence. Simply click on and off for easy replacement and cleaning.

The only true sonic technology



Your Sonicare brush head and handle are the ultimate team. Enhancing our sonic technology, they give you unrivaled cleaning. High-speed vibrations in your toothbrush handle power up to 31,000 brush strokes per minute that whip up your toothpaste into plaque fighting bubbles, driving them deep between your teeth and along the gum line for exceptional cleaning.

Brushing with peace of mind



All Philips Sonicare brush heads are safe and gentle on teeth and gums. During tests, each brush head is put through its paces to ensure durability and an exceptional performance every time you brush.

Specifications

Items included

• Brush heads: 2 i InterCare

Compatibility

- Brush head system: Click-on
- Suitable for these models: 2 Series plaque control, 2 Series plaque defense, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, PowerUp, Protective Clean
- BrushSync mode-pairing

Design and finishing

- Bristle stiffness feel: Regular
- · Color: White
- · Reminder bristles: Blue bristle color fade away
- Size: Standard

Health benefits

- Gum health: Improves gum health in 2 weeks
- Plaque removal: Removes up to 7x more plaque*

Quality and performance

- Replacement: Every 3 months
- Tested: for optimal usage



Issue date 2017-12-01

Version: 2.0.1

© 2017 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

* than a manual toothbrush **BrushSync™ mode-pairing is only compatible with Philips Sonicare BrushSync™ enabled toothbrush handles